

**3M** Science.  
Applied to Life.™

3M™ RelyX™ Unicem and Unicem 2  
Self-Adhesive Resin Cements



Cements:  
Self Adhesive

Long-term  
performer:  
Cement

# All-round reliability.



**3M™ RelyX™  
Unicem**  
Cements


## Clinically proven self-adhesive resin cement

For simplicity, convenience and reliability, 3M™ RelyX™ Unicem Self-Adhesive Resin Cement is your 'go-to' product for a wide range of indications. Thanks to the dual-cure, self-adhesive chemistry, Unicem cements deliver high adhesion, resistance to discolouration, long-term stability and more.

### Ideal for:

- ▶ general cementation
- ▶ oxide ceramic restorations
- ▶ posts and bridges

## Results of *The Dental Advisor* 15 year clinical performance study



Source: *The Dental Advisor* 15-year Clinical Performance Study, March 2017  
Rating system: Excellent +++++

# Formulated to meet the needs of your practice and your patients!

- ▶ Saves time by eliminating etching, priming and bonding steps
- ▶ Virtually no post-op sensitivity
- ▶ Strong
- ▶ Moisture tolerant
- ▶ Endo tips make cementing posts easy and reliable



## 3M™ RelyX™ Unicem 2 Automix Refills

- ▶ Ideally suited for inlays, onlays, crowns and posts
- ▶ Hygienic unidose
- ▶ Consistent mix with any mixing machine
- ▶ 3M™ RelyX™ Unicem 2 Automix Elongation Tip for virtually void-free cementation of posts
- ▶ Working time: 2 minutes



## 3M™ RelyX™ Unicem 2 Clicker™ Dispenser

- ▶ Ideally suited for multiple-unit bridges
- ▶ Hygienic unidose
- ▶ Consistent mix with any mixing machine
- ▶ Working time: 2.5 minutes



## 3M™ RelyX™ Unicem 2 Automix Mix Tips

- ▶ Suitable for all sizes of restorations
- ▶ Choose amount dispensed – 8.5g syringe gives approximately 15 applications
- ▶ Simple and time-saving one-step solution with no need for mixer, activators or appliers
- ▶ Comes with a regular tips for dispensing into crown, and wide tips



## 3M™ RelyX™ Unicem 2 Automix Dispenser

- ▶ For all sizes of restorations
- ▶ Choose amount dispensed – 11g dispensed gives 40 applications
- ▶ Time-saving, easy and hygienic handling
- ▶ Delivers pre-measured doses for consistent and economical mixing
- ▶ No need for mixer, activators, appliers, mixing tips or devices

## Recommended applications ( ++ Highly recommended / + Recommended)

Indication	Metal/ Metal based	Glass ceramics (including Li. Disilicate e.g. e.max®, VITA Mark II)	Oxide ceramics (e.g. 3M™ Lava™ Plus, Brux Zir®, Procera®)	Resin composites (e.g. VITAENAMIC®)
Inlays/Onlays	++	+	++	+
Crowns/Bridges	++	++	++	++*
Endodontic posts	++		++	++
Maryland bridges	+		+	
On implant abutments	++	+	++	++

\*crowns only

## Ordering Information

### Item no. Product information

3M™ RelyX™ Unicem 2 Automix Refills		
56848	Shade Translucent	1 × 8.5g syringes
56846	Shade A2 Universal	10 × regular mixing tips
56847	Shade A3 Opaque	5 × wide mixing tips 5 × endo tips

3M™ RelyX™ Unicem 2 Automix Mix Tips	
56851	30 mixing tips (regular)
56918	15 mixing tips (wide); 15 intra-oral tips
56919	15 mixing tips (wide); 15 endo tips

3M™ RelyX™ Unicem 2 Clicker™ Dispenser Refills		
56874	Shade Translucent	1 × 11g Clicker dispenser
56875	Shade A2 Universal	1 × mix pad (70144)
56876	Shade A3 Opaque	Technique card Instructions for use

3M™ RelyX™ Unicem – Aplicap™ and Maxicap™ capsules	
56828	Trial pack 20 Aplicap capsules (10 × A2 Universal, 6 × Translucent, 4 × A3 Opaque)
56829	Trial kit 20 Aplicap capsules (10 × A2 Universal, 6 × Translucent, 4 × A3 Opaque), Aplicap activator, Aplicap applier
56814	Intro kit 50 Aplicap capsules (24 × A2 Universal, 16 × Translucent, 10 × A3 Opaque), Aplicap activator, Aplicap applier

### Item no. Product information

3M™ RelyX™ Unicem refills	
56815	50 Aplicap capsules (24 × A2 Universal, 16 × Translucent, 10 × A3 Opaque)
56816	50 Aplicap capsules, Translucent, 10 elongation tips
56817	50 Aplicap capsules, A3 Opaque, 10 elongation tips
56818	50 Aplicap capsules, A2 Universal, 10 elongation tips
56819	20 Aplicap capsules, A1
56820	20 Aplicap capsules, White opaque

3M™ RelyX™ Unicem – Maxicap™ capsules	
56830	Intro kit 20 Maxicap capsules (10 × A2 Universal, 6 × Translucent, 4 × A3 Opaque), Maxicap activator, Maxicap applier

3M™ RelyX™ Unicem – Maxicap™ capsules	
56831	20 Maxicap capsules (10 × A2 Universal, 6 × Translucent, 4 × A3 Opaque)
56832	20 Maxicap capsules, Translucent
56833	20 Maxicap capsules, A3 Opaque
56834	20 Maxicap capsules, A2 Universal

3M Oral Care  
3M United Kingdom PLC  
Charnwood Campus  
10 Bakewell Road  
Loughborough  
Leicestershire LE11 5RB  
0845 602 5091

3M Ireland Limited  
The Iveagh Building  
The Park  
Carrickmines  
Dublin 18  
Ireland  
www.3M.co.uk/oralcare

GD025. J394376.  
e.max is a registered trademark of Ivoclar Vivadent.  
Brux Zir is a registered trademark of Glidewell Laboratories.  
Procera is a registered trademark of Nobel Biocare.  
Vita Enamic is a registered trademark of Vita Zahnfabrik.  
3M, Aplicap, Clicker, Lava, Maxicap and RelyX are trademarks of the 3M company.  
© 3M 2018. All rights reserved.

